



Living up to Life




TCS SP8
FSU
Summery Manual
LAS X 3.1

20170703


Contents

I. Start up	2
II. Image acquisition	6
III. Viewer	20
IV. MaxProjection	23
V. Quantify menu	24
VI. Data saving	25
VII. Shut down	30
VIII. Microscope operation	32
IX. HyVolution (Cofocal super resolution method)	36

I. Start up



Main Switches



Indicator lamp of imaging lasers

④ Key switch of imaging laser

① Switch for PC

② Power switch for the Scanner

③ Main power supply of imaging lasers and cooling

- ① Switch on the PC.
- ② Switch on the Scanner.
- ③ Switch on the main power supply of imaging lasers.
- ④ Switch on the key switch of the imaging lasers (yellow indicator lamp will be on)
- ⑤ Switch on the HG lamp.


Note) Do not switch on/off frequently and wait about 5min when restart.



- ⑥ Switch on the microscope.



- ⑦ Log in as "TCS User"
- ⑧ Double click on the "LAS X" on the desktop.




⑨ The following window comes up.

Select modes/switches and click “OK” and wait until the interface opens without operation including microscope.

Configuration: machine.xlhw → Imaging through the microscope.

SimulatorSP8.xlhw → Starting up without connected hardware(Simulator mode). Turn on PC only.

Microscope: DMi8



Select the Scanner
OFF:
FOV scanner; High resolution


ON: Option
Resonant scanner; High speed

⑩ Switch on the lasers

Open the Configuration tab, click the Laser icon.





Activate each laser by single click (Do not double click.)



- Set **Ar** laser power as 20-30%.
- **WLL** laser power will be set as 70% automatically, do not change the power.


⑪ **Change the bit depth**

Click the Configuration menu , and click the Hardware . And select an appropriate number of bit depth.




⑫ **If using Resonant mode; Define the zoom factor**


With Resonant mode, the zoom factor set as bigger when the software opened (zoom factor depends on the 8/12kHz.). Define the zoom factor to minimum before starting acquisition. (Default zoom factor with FOV scanner is 1.0 and minimum is 0.75)




⑬ **if using Resonant mode; set the Live averaging**

In order to get higher image quality with resonant mode under Live scan, set the Live" Line Average during Live Acquisition" in the "Hardware" in the Configuration tab .


The live image will be averaged by setting number for Line average factor. (See II. Image Acquisition)



II. Image Acquisition



Open Acquire Menu



The screenshot shows the 'Acquire' menu in the software. A red arrow points to the 'Acquire' button in the top navigation bar. A callout box explains the 'Wavelength settings' section, which includes selecting lasers, detectors, emission wavelength, and objective lenses. Another callout box points to the '3 buttons for scanning' (Live, Capture Image, Start) at the bottom of the interface. A third callout box points to the 'Viewer' area, which displays a 2x2 grid of fluorescence images (blue, green, red, and merged).

Wavelength settings:
Selecting lasers, detectors, emission wavelength, objective lenses. To simplify the operation, load the defined settings.

3 buttons for scanning

Viewer:
Display area. see IV.Viewer for detail.

Switching Acquisition/Projects tab

Acquisition : Definition of acquire parameters. Pixels, zoom factor, Z stack settings, etc.

Projects : Data container
Acquired image data using by Capture Image button or Start button are

3 buttons for scanning

Live : for preview/adjust button. Data will not be stocked.

Capture Image : for single image acquisition.

Start : for multiple image acquisition.

Wavelength settings

The screenshot shows the 'Acquire' tab of the Leica software interface. The top navigation bar includes 'Configuration', 'Acquire', 'Process', 'Quantify', and 'Analysis'. Below this, there are options for 'Load | Save | Roi' and 'Load/Save single setting: Leica Settings'. A callout box points to the 'Leica Settings' dropdown with the text 'Settings for single color'. To the right, there are checkboxes for 'ROI: OFF', 'Set Background: OFF', and 'Bleachpoint: OFF'. The main area features a color spectrum bar with three vertical sliders. Callouts point to these sliders with the text 'Laser wavelength selection and power'. Below the spectrum bar, there are fields for 'Objective: HC FLUOTAR L 25x/0.95 WATER' and 'Fluo Turret: Scan-BF'. A callout box points to a panel of icons with the text 'Panel box, notch filter, Dye assistant settings'. Below this is a 'Specimen' label and a graph showing emission spectra. Callouts point to the graph with the text 'detectors (PMT or HyD)' and 'Detection setting (Prism-slit)'. Below the graph, there are settings for 'HyD 1' (ON, Gain 100.0), 'PMT 2' (OFF, Gain 800.0), and 'HyD 3' (OFF, Gain 100.0). At the bottom, there is a 'TLD' section with 'PMT Trans: OFF' and 'Gain [V]: 160.4'. A callout box points to this section with the text 'Transmission detector (Bright field)'.

① **Select objective**


Select objective and observe specimen, find the position and the focus.



① **Load the fluorescence setting For Single staining observation**


Select the light path setting from the pull-down menu of “Leica Settings” of “Load/Save single settings”. *For the multi staining, go to ③.

DAPI: UV ex. Blue fluorescence observation
FITC: Blue ex. Green fluorescence observation (Alexa488, Cy2, GFP etc)
TRITC: Green ex. Red fluorescence observation (Alexa555,568, Cy3, DsRed etc)
Cy5: Red ex. Far Red fluorescence observation (Alexa633,647, TOTO3, TOPRO3 etc)
Transmission: bright field




② **Load the fluorescence setting For multiple staining observation**


③ – 1 Activate “SEQ”button in the “Acquisition mode” panel. “Sequential Scan” panel will open.



③ – 2 Click the “Load”button on the “Sequential scan” panel and select the sequential setting for the multiple staining.



③—3 Several “Seq.” buttons are displayed depend on the number of fluorophores.
 Activate each “Seq”button and define for the each staining (laser power, gain).
 Refer ④ and after how to define these parameters.




*** Sequential mode**

between lines		Each excitation changes in line by line. Suitable for live imaging since it is almost simultaneous scanning, time lag is only a line period. Necessary to have same number of detectors with dyes.
between frames		Image acquisition executes each frame in different excitation. It means there is a time lag of a frame period.
between stacks		Each excitation changes in stack by stack. This method achieves the fastest completion time even the time lag between the dyes becomes bigger

Loading the settings from the acquired/saved data

By using the Apply button on the Experiments tab, the settings of the acquired/saved data can be loading (reproducible). Activate the raw image data on the Projects tab and click Apply. All the parameters of wavelength(Excitation, emission), detectors will be reproduced.



④ Image acquisition

Click the Live button  at the lower left on the monitor.

Scan starts and align the parameters as follows.

Notes: Live button expose lasers your specimens. So it should be stopped often in order to avoid the bleaching.

■ Defining the laser output

Define the laser output so that it may be optimized to the specimen by using slider or inputting the number on the laser panel.

Switching WLL and Ar laser, 405nm laser.


Displays switch by clicking the vertical buttons "Switch to Whitelight", "Switch to Conventional"



Defining laser power by using slider or inputting numbers.



* Activate corresponding Seq button whenever define each laser power.




First activate corresponding Seq button and change the laser power.

Defining the focus, and the detector gains

Adjust focus and the detector gains and offsets by using each dial of the control panel (option) as below.

For the multistaining, the “Smart Gain” and the “Smart Offset” dials can adjust each staining by activating channel image respectively on the monitor.

- default settings of gains as a standard → PMT:800V HyD:100%




- point

The less Gain achieves the higher S/N. If higher gain is used, use “Average” function in order to reduce noises.

	PMT (conventional)	HyD (Hybrid detector)
Gain (the voltage of detector) Default Gain as a standard	0-1250V 800V	10-500% 100%
Offset (background setting)	Possible	-
Photoncounting mode	-	possible

For multicolor imaging, with Panel Box, you can change each channel Gain/Offset by using Smart Gain/Smart Offset dial, just click certain (singlecolor) image on the monitor,





Quick Look Up Table (Q LUT) button

For easy tuning of signal intensity, click Q LUT button and the look up table (image display color) changes to “Grow Over/Under”.

Image display color changes like follows by clicking this button.

Original (pseudo color like green, red, etc.)

→Grow Over/Under (for checking)

→Gray (black white)

→Original (pseudo color like green, red, etc.)

Image can be taken with 8/12/16bit. For instance 8bit; it has 256 intensity levels (grey scale/gradation), and “Grow Over/Under” color shows this 256 level in different colors like follows.



>255	blue (saturation = overflow)
255	blue
·	white
·	brown
·	black
0	green
0 <	green

Blue color means the signal saturation. Readjust laser power and/or gain to avoid saturation.

Green color means the darkest signal. Too low offset makes enlarged green area and some weak signals may be failed to be detected.

Green:
Background
intensity = 0.


Blue:
saturation.



⑤ Other parameters (format, scan speed, zoom, etc)

If necessary, Format (pixel numbers of the acquired image), Speed (scan speed), Zoom factor (optical zoom) can be also changed. The slower the scan speed, the smoother the image quality.

Please note, the bleaching effects becomes bigger with the higher the zoom factor.



Format: Pixel numbers of the image

Speed: scan speed

Zoom: minimum is 0.75.
Zoom In button enables you to define the region of interest on the image.
Control panel also can be used.

Rotation(option)
Scanning rotation

Zoom
Zoom factor

Average: averaging for reducing noise.

Rotation(option) : scanning rotation

⑥ Average function

“Average” enables to get higher signal to noise (image quality) ratio of the image.

- Set higher number of average for higher detector Gain.
- Higher number of average does not necessary for lower Gain.

⑦ Acquisition of the Single image


Click “Capture Image button”  for single slice image.


Image data will be kept in the Project tab as the name of “Image001” temporarily.
The number of Average will be activated by these buttons during the acquisition.

⑧ Z stack acquisition

Define start position with “Begin” button and end position with “End” button while scanning by “Live” button. Focus can be moved by Z position dial of the Control panel.



Z position
Adjust z position/focus



⑧-1
Define Begin/End position

Delete defined Begin/End position

Thickness by defined Begin/End

⑧-2
Define number of frames or step size.

Guideline for z-step size: About the twice of the number of “System Optimized”

Activate either of them

Z-Stack : 192.82 μ m | 558 Steps

Begin [μ m] : 95.15

End [μ m] : -97.68

Z Position [μ m] : -6.00

Z Size [μ m] : 192.83

Re-Center

Stack Direction (Z) : \uparrow \downarrow

Nr. of Steps : 558

Z-Step Size : 0.35

System Optimized

Z-Compensation : AOTF/EOM & Detector

Galvo Flow : OFF

Travel Range [μ m] : 500

⑨ Acquisition of the Zstack

Data acquisition starts by clicking following button.

XYZ (or for several continuous image acquisition:  click “Start” button.


Image data will be kept in the Project tab as a name of “Series001” temporarily.

Note: keep away from the microscope and the anti-vibration table while scanning in order to avoid the vibration noise to the images.

Parameters of the acquired image data

Right mouse click the data on the Experiments tab and select “Properties”.

“Experiment Data” will open and able to confirm the parameters of the image data.



The image shows two parts of the software interface. On the left, a context menu is open over 'Image001' in the 'Projects' list. The 'Properties' option is highlighted with a dashed red box, and a grey arrow points from it to the right. On the right, the 'Properties' dialog box is displayed, showing detailed acquisition parameters for 'Image: Series004'.

Properties Dialog Box Data:

Description: [Apply Text]

Image: Series004
 Size: 6.29 MB
 File Locator:
 Start Time: 20120704 17:40:01.125
 End Time: 20120704 17:40:01.163
 Total Exposures: 24 (1 channels, 24 frames)

Dimensions			
Dimension	Logical Size	Physical Length	Physical Origin
X	512	250.00 µm	0.00 µm
Y	512	250.00 µm	0.00 µm
Z	24	-11.23 µm	-0.02 µm


Channels			
LUT	Resolution	Min	Max
Blue	8	0.000000e+000	2.550000e+002

Time Stamps			
Frame (Ubuu:All)	Relative Time (s)	Absolute Time (h:m:s.ms)	
1	0	17:40:01.125	:
24	0.058	17:40:01.163	:

Confocal Settings	
Name	Value
Scan Mode	xyz
Scan Direction X	Unidirectional
Scan Speed	400 Hz

HyD (Hybrid Detector *option)

HyD has higher sensitivity, higher S/N, and extremely low noises. Photoncounting mode can be selected.



Caution:

- Please separate mobile items (cellar phone and other items which send electric wave) from HyD about 1m. If they are close to locate each other, HyD would detect noises from them.
- If HyD detects extremely high intensity, its shutter closes automatically. It means too much emission fluorescence by high excitation power. Redefine lower laser power.

The difference from PMT

- HyD Gain is variable from 10 to 500% whose PMT is variable from 0V to 1250V.
- NO Offset function for HyD.
- There are three kinds of modes with HyD.

Detection mode for HyD

- **Standard:** Normal mode as PMT.
- **Counting:** Please refer next page.
- **BrightR:** Mode for imaging with wider dynamic range a little bit.



Selection of HyD modes

Counting mode (Photon Counting mode)

This is a mode to count and accumulate photons in each pixel. In this case, the intensity means the number of photons. Detector noises are almost zero and the S/N and the quantitativity improves quite a lot.


Theoretically it is not possible to use Gain with Counting mode, the image shows less bright than Standard mode, to get brighter image with Counting mode, please refer following settings.

<points of using Counting mode>

- Gain cannot be used in theoretically. Only the number of photons are shown as a image.
- brightness adjustment can be done with following method.


① take slower scanning speed. (only with FOV scanner)

slower scanning speed can get more photons per pixel which means brighter image.



② take Accumulation

For example, 2 times Accumulation made about twice as bright image.




⑩ **timelapse**

⑩—1. Select mode

Select “xyt” or “xyzt” from the Acquisition Mode.


xyt : XY timelapse imaging without Z stack

xyzt : XYZ timelapse imaging with Z stack.



⑩—2. Following panel will open according to the selected mode.


Define the “Duration”(whole acquisition time) and “Time interval”



Frames/Stacks will show depend on the mode(xyλ/xyzt). If the xyzt mode is selected, this column will be Stacks and Z settings should be defined on the Z stack panel.

⑩—3. Timelapse acquisition



Define other settings (image adjustment, zoom, Z stack, etc.) and click “Start” button for scanning.




III. Viewer

Viewer: Window which is shown the acquired images. There are function buttons around the images.










The colors of the image can be changed by clicking this LUT bar.









Tiled/full display can be changed by double clicking the image.




Icons of the viewer

	<p>Select tool</p>
	<p>Delete selected/activated item</p>
	<p>Insert scale.</p> <p>Activate the scale by using  button and right mouse click, select "Properties" then define the length, angle, etc.</p> 
 151.95 %  100 %	<p>Display scale.</p> <p>On/Off this button and the image changes 100% display or shrink or zoom depend on the image resolution and the display resolution.</p>
	<p>Quick Look Up Table</p> <p>The color of the image changes by clicking this button each time, Grow Over/Under→Gray→Original.</p>  <p>Original color</p> <p>Grow over/under Blue : indicates saturation</p> <p>Gray</p>


	<p>Overlay Display or not display “overlay image” of the multiple channels by activating/not activating this button.</p>
	<p>Maximum projection Create Maximum projection image from the XYZ stack data temporarily. The maximum intensity in the z stack for each pixel will be shown as a single 2D image.</p> 
	<p>Orthogonal Sectioning Show XZ and YZ section images from the XYZ image data.</p> 
	<p>3D 3D Viewer launches with this button. 3D view will be shown from the XYZ data.</p>

IV. MaxProjection

Create Maximum Projection image data from the Z stack image.







Open the "Process" menu and select the "Projection" in the Process Tools. Activate the Z stack data on the Experiments tab and click "Apply" button. Maximum Projection image data will create/add in the Experiment tab.



V. Quantify menu

Functions to measure lengths, intensities, areas, histograms, ratio, etc.
There are three mode as follows.



	Line	Lengths and intensities and other statistical data are shown by drawing line on the images. This function is supposed to use for single image.
	Stack	Intensities, areas and other statistical data of defined ROIs (Region of Interest) are shown. This function is supposed to use for XYZ or timelapse data.
	Histogram	Histogram and the other statistical data are shown by defined ROIs (Region of interest). This function is supposed to use for single image.

Result export by “Report” button


Result sheet can be export by “Report” button.










Saved data by this button will be .xml, jpg, csv files. Graph data will be imported by Excel, as well.

VI. Data saving

All image data in the single file will be saved as a .lif file.




Project Tool bar


-  **New:** Create new Project file.
-  **Open:** Open saved data.
-  **Save:** Save all data on the Projects tab.
-  **Apply:** Restore parameters of which activated data. Selection of applying parameters can be defined from following mene. Configuration > Instrument Parameter Settings (IPS)
-  **Browse:** Browsing saved data.
-  **Shows:** Show the first image data on the Projects tab.

① **Right mouse click the data and select "Save As".**


Image data in the same file will be saved as a single file as explained above.

All image data in the same file will be saved as a single .lif file.

To rename each data (Series001,Image001, etc), right mouse click and select "Rename" or just double clicking the data.




② **Select the path to save the data.**



③ **Overwrite saving (refresh saving)**

Adding/editing data after saving .lif file, data should be overwritten by executing “Save Project”



④ **Acquiring new data in the new file**

To create new file, click “New” button, new data will be acquired in the activated new file.


New button

Created new file
 Activate this file and acquire new data, the data will be input in the new file.
 * data will be able to copy/move by click and drop/copy and paste to the other file.


Note)
If the new data are acquired without activating new file, it will be acquired in the (another) activated file.

⑤ Exporting file as “.Tiff”, “.Jpg” format

Right mouse click the data and select "Export" and select "As Tiff..." or "As JPEG..".



Then following dialog will open.
Select the path and click save.




LAS X Download site

1. Access to following site.

<http://www.leica-microsystems.com/products/microscope-software/software-for-life-science-research/las-x/>


2. There is a link at the most end of the page, click to download.





The screenshot shows the Leica Microsystems website's 'Downloads' section. At the top, there is a navigation bar with links for HOME, PRODUCTS, APPLICATIONS, NEWS & MEDIA, COMPANY, CONTACT, SEARCH, and QUOTATION CART: 0. The main content area is dark grey with a 'Downloads' heading. A dashed white box highlights three download links: 'Download free LAS X Core Offline version 1.1.0 for Windows7 64Bit', 'Download free LAS AF Lite version for Windows XP, Windows7 32Bit, Windows7 64Bit', and 'Download Leica LAS X Brochure (English)'. Below this is an 'Available on the App Store' badge and links for 'Download the LAS X Mobile Connection app free of charge' and 'Download Leica SDK Hardware Configurator version 2014.3'. A large white arrow points from the highlighted links to a second 'Downloads' section below. This second section lists three items: '1 • Download free LAS X Core Offline version 3.0.1 for Windows7 64Bit', '2 • Download free LAS AF Lite version for Windows XP, Windows7 32Bit, Windows7 64Bit', and '• Download Leica LAS X Brochure (English)'. A callout box on the right of the second section contains the text '1 For SP8'.

VII. Shut down

1. Switch off the lasers on the software.

Open the Configuration tab  and click “Laser Config”.



- Ar laser: Decrease the laser power down to 0% and switch off.
- WLL laser: Switch off while keeping laser power 70%

2. Lens cleaning

Using lens paper/ cotton-tipped stick and lens cleaner, clean up liquid immersion lens.

Note) Do not use Kim wipe or other papers for lens cleaning.

3. Close LAS X and Windows.

Save all data and close the LAS X and shut down the Windows.

4. Switch off the mercury lamp .

Note; Wait about 5min. for restart.



5. Switch off the Microscope.



6. Switch off main switches

- ① If the Windows is already finished, switch off the PC.
- ② Switch off the scanner
- ③ (If Ar laser is equipped) after the cooling fan was stopped (about 5min after deactivate Ar laser on the software), switch off the key switch.
- ④ Switch off the laser power supply.



③ Key switch of imaging lasers.

① Switch for PC


② Power switch for the Scanner

④ Main power supply of imaging lasers and cooling fan of Ar.laser

VIII. Microscope operation


1. Focusing and XY Stage control



Focus the image by turning the focus dials on the left and right sides of the stand. Alternatively, rear rotatory knob on the Leica Smart Move control element can also be used. XY stage can be controlled by Leica Smart Move control element.



2. Changing objective nose pieces

2-1 Changing objective by using the Touch Screen




- ① Use the  tab to switch to the Magnification menu. Currently selected objective is displayed in red.
- ② To switch between objectives, press the corresponding objective key from Objective Nosepiece Mode. Immersion objectives are marked with drop symbol  63. When you change the objective between IMM and Dry, you need to touch the objective key

two times. The objective key will blink by the first touch and the objective will be changed by the second touch.

Note) If the immersion liquid is attached to dry lens, clean the dry lens certainly. Otherwise it causes the decrease of image quality.

2-2 Changing objective on the software.

Click 'Objective' in the software and select the appropriate lens from the pull-down menu.



※ The message is shown when changing between dry lens and immersion lens.

“Immersion of current objective(OIL/IMM/DRY) is different from selected objective(OIL/IMM/DRY).
Do you want to turn turret automatically? Yes No ”

→ **Crick Yes** to change the objective lens automatically

In changing from Immersion lens to Dry lens, wipe off the Immersion from the specimen to prevent Immersion forming on the objective lens.

※ Immersion (IMM) Objectives

For immersion objectives use the appropriate immersion medium.

OIL: Use optical immersion oil from Leica only

W: Water immersion


Gly: glycerin


IMM: Universal objective for water, glycerin, Oil immersion. (Require setting of the appropriate correction circle to use each IMM)

Note) If the immersion liquid is attached to dry lens, clean the dry lens certainly. Otherwise it causes the decrease of image quality.

3. Microscope observation method (by your eyes.)

3-1. Transmitted method(Including DIC)



① On the Touch Screen use the  tab to configure the contrast method.

② Select the contrast method from Transmitted method list.



BF: Brightfield Transmitted Light


DIC: Differential Interference Contrast

③ Touch the TL-Shutter key on the Touch Screen to open or close the Shutter (active illumination axis on or off). Alternatively, Shutter control button on the left side of the stand can


also be used.


④ The brightness is adjusted by the knob on the left side of the stand.




⑤ On the Touch Screen use the  tab to confirm the current intensity.

3—2. Operating the fluorescence





- ① On the Touch Screen use the  tab to configure the contrast method.
- ② Select FLUO
- ③ Select the desired cube using the corresponding key from FLUO-Filtercubes.
※ Set of Filter cubes are depending on the system. Please refer the Fig1.1 about the spec of each FLUO filter cubes.
- ④ Touch the IL-Shutter key on the Touch Screen to open or close the Shutter (active illumination axis on or off). Alternatively, Shutter control button on the left side of the stand can also be used.

⑤ The brightness is adjusted by the knob on the left side of the stand.



Filter	Typical fluorophore	Excitation	Excitation filter	Dichroic	Fluorescent filter
DAPI	DAPI	UV	BP350/50	400	BP 460/50
FITC	FITC (Band Pass)	Blue	BP480/40	505	BP527/30
FITC LP	FITC (Long Pass)	Blue	BP470/40	510	LP515
RHOD	Rhodamine	Green	BP515-560	580	LP590



- ⑥ On the Touch Screen use the  tab to confirm the current FIM (Fluorescence Intensity Manager). The brightness is adjusted in 5 defined increments.

Leica Microsystems K.K.

Tokyo

1-29-9 Takadanobaba, Tokyo, 169-0075 Japan

Phone: +81-3-6758-5630

Fax: +81-3-5155-4333

Osaka

10F 5-4-9 Toyosaki, Kitaku, Osaka City, Osaka 531-0072 Japan

Phone: +81-6-6374-9771

Fax: +81-6-6374-9772

Nagoya

2F 2-15-20 Nishiki, Nakaku, Nagoya City, Aichi, 460-0003 Japan

Phone: +81-52-222-3939

Fax: +81-52-222-3784

Fukuoka

12F, 8-30 Tenyamachi, Hakataku, Fukuoka City, Fukuoka 812-0025 Japan

Phone: +81-92-282-9771

Fax: +81-92-282-9772



ライカ マイクロシステムズ 株式会社

本社 〒169-0075 東京都新宿区高田馬場1-29-9
Tel.03-6758-5640 Fax.03-5155-4336

大阪セールスオフィス 〒531-0072 大阪府大阪市北区豊崎5-4-9 商業第2ビル
Tel.06-6374-9771 Fax.06-6374-9772

名古屋セールスオフィス 〒460-0003 愛知県名古屋市中区錦2-15-20 三永伏見ビル2F
Tel.052-222-3939 Fax.052-222-3784

福岡セールスオフィス 〒812-0025 福岡県福岡市博多区店屋町8-30 博多フコク生命ビル12F
Tel.092-282-9771 Fax.092-282-9772

●<http://www.leica-microsystems.co.jp> E-mail: lmc@leica-microsystems.co.jp

※この製品のデザインおよび仕様は改良などのために予告なく変更する場合があります